

FASTER ORTHODONTIC TREATMENT

AcceleDent® Aura is a simple-to-use, hands-free device designed to accelerate tooth movement - with only 20 minutes daily use.



Fast

- ›› AcceleDent Aura advances the science of orthodontic care with innovative SoftPulse Technology™.
- ›› Safe acceleration of the bone remodeling process complements conventional orthodontic treatment.


AcceleDent Aura delivers compelling clinical results

- ›› Clinically proven to accelerate tooth movement up to 50%.¹
- ›› SoftPulse Technology is based on demonstrated scientific methods in orthopedic and orthodontic bone remodeling.^{2,3}

Safe

Demonstrated safe and reliable in U.S. clinical trials

- ›› Investigators reported no risk of increased root resorption and no risk in increased loss of posterior anchorage.^{1,3}
- ›› AcceleDent Aura is a Class II medical device cleared by the FDA.



SoftPulse Technology™ generates precisely calibrated micropulses to accelerate remodeling activity during orthodontic treatment.

AcceleDent®
aura

"AcceleDent gives our patients an attractive option to have the smile they desire in a time frame they will accept."

– Robert Rudman, DDS, MS



Gentle

AcceleDent may improve the orthodontic experience for patients.

» SoftPulse Technology micropulses are gentle and exert up to 8x less force than a power toothbrush.




■ AcceleDent Aura

vs.



■ Oral-B® Sonic Complete power toothbrush



The small, lightweight **Activator** generates gentle micropulses to safely accelerate bone remodeling. It includes a USB interface which can be plugged directly into a computer to review patient usage history via the FastTrac Report.



The **Mouthpiece** provides a comfortable fit and snaps easily on and off for transport and cleaning.



The **USB Extension Cable** with **Power Adapter** allows for easy, convenient charging.



The new **Travel Case** keeps the components safe, clean and secure.



"My patients and I are loving AcceleDent. With the shortened treatment time I am convinced that our patients will be better focused throughout treatment, seeing results sooner."

– David J. Henderson, DDS

AcceleDent
aura

References:

1. Data on file.
2. Rubin C, Judex S, Qin YX. Low-level mechanical signals and their potential as a non-pharmacological intervention for osteoporosis. *Age Ageing*. 2006 Sep; 35 suppl 2:ii32-ii36.
3. Kau CH. A radiographic analysis of tooth morphology following the use of a novel cyclical force device in orthodontics. *Head Face Med*. 2011 Aug 9; 7-14.

AcceleDent is an orthodontic accessory intended for use during orthodontic treatment. Refer to the Directions for Use with AcceleDent Aura or acceleddent.com for warnings and precautions prior to using this product. Federal (U.S.) law restricts this device to sale by or on the order of a licensed doctor. Rx use only for a single patient.

Made in the USA

©2013 OrthoAccel® Technologies, Inc.

AcceleDent® is a registered trademark of OrthoAccel® Technologies, Inc.

Oral-B® is a registered trademark of P&G.



L60000 Rev.2, JUNE 2013

